VIPERSLIDE® Lubricant

Instructions for use

The following are trademarks of Cardiovascular Systems, Inc.:

CSI®, Cardiovascular Systems®, VIPERWIRE Advance®, VIPERSLIDE®

Explanation of symbols on the package and package label

	Conformité Européenne (European Conformity)
	EU Authorized Representative
	Model number
	Caution: Consult IFU www.csi360.com (USA)
	Sterilized with steam
	Lot number
	Do not reuse
	Do not re-sterilize
	Manufacturer
	Use by
	Temperature limits
	Fragile
	This side up
	Keep dry
	Keep out of direct sunlight

1. Description

The Cardiovascular Systems, Inc. (CSI) VIPERSLIDE Lubricant is designed to reduce the friction between the flexible drive shaft of the CSI Orbital Atherectomy Device and the ViperWire Advance Guide Wire. The VIPERSLIDE lubricant is an emulsion composed of soybean oil, egg yolk phospholipids, glycerin, and water for injection. An addition of sodium hydroxide adjusts the pH resulting in a final pH of 8.0 with a pH range of 6.0 to 8.9. See the table below for a list of ingredients

and amounts per ingredient. The VIPERSLIDE lubricant is provided in a sterile, non-pyrogenic 100 mL bag.

Ingredient	Amount
Soybean Oil	10 %
Egg Yolk Phospholipids	1.2 %
Glycerin	2.25 %
Sodium Hydroxide (pH range is 6.0 to 8.9)	QS*
Water for Injection	QS*

QS* = Quantity Sufficient

2. Contents of package

Each package contains the following:

- 10 - 100 mL **single-use** bags of VIPERSLIDE Lubricant

3. Indications

VIPERSLIDE Lubricant is indicated for use with CSI Orbital Atherectomy Systems.

4. Contraindications

VIPERSLIDE Lubricant is contraindicated in patients with hypersensitivity to egg-, soya-, or peanut protein or to any of the active substances or excipients of the lubricant components listed under Section 1 above.

5. Warnings

- Do **not** store or use partially-used bags of lubricant as the lubricant is for **single-use** only. If VIPERSLIDE is re-used, it may lead to serious infection and/or patient harm including death. Dispose of all single-use lubricants according to local environment requirements.

- Do **not** expose or store the lubricant package to below freezing temperatures. If the lubricant package is stored or exposed to below freezing temperatures, the lubricant will destabilize. Discard the bag of lubricant and obtain a new bag of lubricant.

- Do **not** use the lubricant if the bag or seal is torn, is punctured, or has leaked. Discard the bag of lubricant, as sterility may be compromised, and obtain a new bag of lubricant.

- Do **not** use the lubricant in which there appears to be an oiling out (separation of oil and water) on the surface of the emulsion. Visually inspect the lubricant for the presence of particulate matter and discoloration prior to use. If oiling out, particulate matter, or discoloration exist within the emulsion, discard the bag of lubricant and obtain a new bag of lubricant.
- Do **not** use the lubricant if the integrity indicator is black. If the integrity indicator is black, discard the bag of lubricant and obtain a new bag of lubricant as the overpouch is damaged.

6. Precautions

- The prime stabilizers of emulsions are excessive acidity (low pH) and inappropriate electrolyte content. Careful consideration should be given to additions of divalent cations (Ca++

and Mg++) which have been shown to cause emulsion instability.


- Consult a pharmacist prior to mixing additives to the VIPERSLIDE/saline admixture. Do not add additives directly to the VIPERSLIDE lubricant.
- The drug Papaverine should not be added to the VIPERSLIDE/saline admixture.
- Prepare the admixture using strict aseptic techniques as the lubricant is a growth medium for microorganisms.
- Begin using the admixture within one hour after admixture preparation. The admixture is considered an immediate-use solution with a Beyond-Use-Date (BUD) of one hour. If administration of the admixture has not begun within one hour, discard the admixture and prepare a new bag of admixture.

7. Storage and handling

- Do **not** re-sterilize or reuse the lubricant bag.
- Use the lubricant prior to the expiration date on the shelf carton label and the individual lubricant bag label.
- Do **not** freeze the bags of lubricant. Discard the bags of lubricant if the bags of lubricant have been frozen.
- Do **not** store bags of lubricant above 25 °C (77 °F).

8. Opening the package

1. Inspect the integrity indicator in the overpouch. If the integrity indicator is black, discard the bag of lubricant as the overpouch is damaged (see callout A below).
2. Tear downward from the notch at the top of the overpouch and remove the lubricant bag (see callout B below).
3. Dispose of the integrity indicator and the oxygen absorber packet (see callout C below).

A Integrity indicator
B Overpouch notch
C Oxygen absorber packet

9. Directions for use

9.1. Preparing admixture with additives

Caution: Do not add additives directly to the VIPERSLIDE lubricant.

Follow the proper mixing sequence below to minimize pH-related problems by ensuring that acidic injections are not mixed directly with the VIPERSLIDE lubricant:

1. Inject additives to the saline bag.
2. Thoroughly mix the additives and saline, with a gentle agitation, to ensure there are no localized concentrations of additives.
3. To complete the admixture, see Section 9.2 below.

9.2. Completing the admixture

Caution: Prepare the admixture using strict aseptic technique as the lubricant is a growth medium for microorganisms.

Note: Refer to CSI's orbital atherectomy device or system instructions for use for preparing for and performing atherectomy procedures.

1. While using aseptic technique, use a syringe to draw 20 mL of lubricant and inject the lubricant into the saline bag. This mixture will be used throughout the atherectomy procedure.
2. Mix at a ratio of 20 mL of lubricant per 1000 mL of saline. Discard any unused portion of the lubricant.
3. Thoroughly mix the lubricant and saline solution by inverting the saline bag three to five times or until the mixture is uniformly opaque. Do not shake the saline bag to mix the lubricant and saline. Use care to mix the lubricant and saline solution in order to avoid forming air bubbles within the solution.
4. Begin using the admixture within one hour after admixture preparation. The admixture is considered an immediate-use solution with a Beyond-Use-Date (BUD) of one hour. If administration of the admixture has not begun within one hour, discard the admixture and prepare a new bag of admixture.
5. Refer to CSI's orbital atherectomy device or system instructions for use for instructions on connecting the bag of admixture to the CSI Saline Infusion Pump or Controller.

10. Disclaimer of warranty

Cardiovascular Systems, Inc. (CSI) disclaims all warranties whether expressed or implied, written or oral, including but not limited to any warranties of merchantability of fitness for a particular purpose. CSI does not warrant either for a good effect or against any ill effect following use. CSI shall not be liable for any direct, indirect, incidental or consequential loss, damage, or expense arising from or related to the use of this device.

CSI may, at its sole discretion, replace any device that is determined to have been out of specification at the time of shipment.

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EC REP **ICON (LR) Limited**