**Diamondback 360 Coronary Orbital Atherectomy System**

Instructions for use

Including the *Diamondback 360* Coronary Orbital Atherectomy Device with *GlideAssist*, Saline Pump, *ViperWire Advance* Coronary Guide Wire, and *ViperWire Advance* with Flex Tip Coronary Guide Wire

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

The following are trademarks belonging to Cardiovascular Systems, Inc.: *CSI*, *Diamondback 360*, *ViperSlide*, *ViperWire Advance*, and *GlideAssist*

**Explanation of symbols on package labels**  
Refer to the package labels to see which symbols apply to specific products.

- Lot number
- Model number
- Consult IFU
- Caution: Consult IFU [www.csi360.com](http://www.csi360.com) (USA) (symbol appears in blue when placed on the device itself)
- Caution: Consult IFU [www.csi360.com](http://www.csi360.com) (USA)
- Do not use if the product sterile barrier system or its packaging is compromised
- Non sterile
- Do not reuse
- Do Not Resterilize
- Manufacturer
- Tip diameter
- Shaft diameter
- Length
- Quantity per box (this example = 5 per box)
- Use by
- Maximum outer diameter
- Crown style
- Nose length
- Serial number
- Sterilized using ethylene oxide gas

- Contains Phthalates
- Contains Phthalates
- Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

**Explanation of symbols on the Saline Pump**

- Low saline red LED indicator
- Start button and green LED pump on indicator
- Pump status yellow LED indicator
- Prime button
- Ground
- Type CF Applied Part
- Slow Blow Type T Fuse
- Finger Caution

**Explanation of Symbols and Graphics on the Coronary Orbital Atherectomy Device**

- Low speed button
- High speed button
- Prime button
- Non-continuous use; spin cycles of 30 seconds on, 30 seconds off with maximum spin time of ≤ 5 minutes

**Table of Contents**

- 1. System description .....1
- 2. Indications for use .....2
- 3. Contraindications .....2
- 4. Warnings .....2
- 5. Precautions .....2
- 6. OAS component storage and handling .....2
- 7. Adverse events .....2
- 8. Clinical trial summary .....2
- 9. Equipment, set up, and test .....2
- 10. OAS Directions Use .....3
- 11. Specifications .....4
- 12. OAS Pump Declaration of Conformity .....4
- 13. EMC Declaration .....4
- Appendices**
- A. OAS troubleshooting .....6
- B. Maximum orbit and result lumen diameter .....7
- C. Orbit performance .....7

**1. System description**

The Cardiovascular Systems, Inc. (CSI) *Diamondback 360* Coronary Orbital Atherectomy System (OAS) is a catheter-based system designed for facilitating stent delivery in patients with coronary artery lesions. The OAS consists of the hand-held *CSI Diamondback 360* Coronary Orbital Atherectomy Device (OAD) with *GlideAssist*, the *CSI Saline Pump* (OAS pump), the *CSI ViperWire Advance* or *ViperWire Advance* with Flex Tip Coronary Guide Wire (*ViperWire* guide wire), and the *CSI ViperSlide* Lubricant. The OAS reduces coronary plaque on the vessel wall by using an orbiting, diamond-coated crown, within coronary arteries, in order to facilitate stent delivery.

- The OAS consists of the following:
- *Diamondback 360* Coronary Orbital Atherectomy Device Model DBEC-125
  - Orbital Atherectomy System (OAS) Pump Model SIP-3000
  - *ViperWire Advance* Coronary Guide Wire Model GWC-12325LG-FLP or *ViperWire Advance* with Flex Tip Coronary Guide Wire Model GWC-12325LG-FT
  - *ViperSlide* Lubricant Model VPR-SLD2

**Caution:** The OAD and *ViperWire* Guide Wire are for use with *Diamondback 360* coronary OAS components only.

**1.1 Orbital Atherectomy Device (OAD)**

The orbital atherectomy device (OAD) is a hand-held, over-the-wire device that includes a sheath-covered drive shaft and a diamond-coated crown (Figure 1). The diamond coating on the crown provides an abrasive surface with which to reduce coronary plaque within coronary arteries. The OAD is designed to track and spin only over the *CSI ViperWire Advance* or *ViperWire Advance* with Flex Tip Coronary Guide Wire. The *GlideAssist* feature facilitates advancing and retracting the OAD crown over the guide wire. **Do not use any other guide wire with the OAD.** See Table 1 for crown size information. See Appendices B and C for orbit performance.

**Table 1. Crown size**


Model Number	DBEC-125
Crown Size (mm)	1.25
Nose Length* (mm)	5
Maximum Lumen (mm)	1.84
OAD maximum catheter outer diameter (mm)	1.45
Minimum guide catheter diameter	6 French with an internal diameter of at least 0.066inches (1.68mm)

\* Nose length is the length of the drive shaft from the crown to the distal tip of the shaft.

- OAD components:
- crown
  - crown advancer knob
  - drive shaft
  - sheath (covering the drive shaft up to the crown)
  - electrical power cord
  - saline tubing for connecting the OAS pump to the OAD

- OAD features:
- On/Off button that allows the user to control when crown rotation starts and stops
  - two (2) speed control buttons on the OAD that allow the user to select the crown rotation speed
  - *GlideAssist* to facilitate advancing and retracting the OAD crown over the *ViperWire* guide wire.
  - crown advancement measurement indicators
  - manual guide wire brake that allows the user to restrict both the rotational and axial movement of the *ViperWire* guide wire

- an eccentrically-mounted, diamond-coated crown that provides an abrasive surface with which to reduce coronary plaque on the vessel wall
- a green polymer strain relief on the nose of the OAD that prevents kinks where the sheath and the drive shaft attach to the handle
- Prime button to momentarily increase saline/lubricant flow

**Figure 1. OAD**

- A. Nose length
- B. Crown diameter
- C. Strain relief
- D. Lockable crown advancer knob
- E. On/Off button
- F. Low and high speed buttons and indicators
- G. Guide wire brake lever
- H. Electrical power cord
- I. OAD saline port connector
- J. Saline bag spike
- K. Saline tubing positioners
- L. Saline tubing
- M. Injection port
- N. Saline line port connector
- O. Sheath
- P. Drive Shaft
- Q. Prime Button


**1.2 OAD package contents**


The OAD and accessories are supplied sterile and are for single-use only. Each package contains:

- OAD Model DBEC-125
- Saline tubing (connects the OAD to the OAS pump)

**1.3 OAS Pump**

The OAS pump provides the saline pumping mechanism and power to the OAD. The small, reusable, and portable OAS pump attaches to a standard five-wheel rolling intravenous (IV) pole (Figure 2) or table mount pole. The OAS pump includes a built-in, audible 25 second spin time notification, OAS power and priming buttons, and status indicators.

**Figure 2.** OAS Pump (Front and Rear)

- A. IV pole screw clamp
- B. IV pole or table mount (not included)
- C. Low saline level sensor and connector cord
- D. Control panel
- E. OAS pump door
- F. OAD connection
- G. Master Power Switch
- H. AC Power Cord Input
- I. Saline Sensor Jack

#### 1.4 OAS Pump package contents

- The OAS pump and accessories are supplied non-sterile. Each package contains:
- OAS pump Model SIP-3000 with attached IV pole screw clamp
  - Power cord
  - Low saline level sensor and connector cord

#### 1.5 ViperWire Advance and ViperWire Advance with Flex Tip Coronary Guide Wire

The *ViperWire* guide wire is a smooth wire, with a silicone coating, and a radiopaque distal spring tip (Figure 3). The *ViperWire* guide wire allows for proper positioning of the OAD crown within coronary arteries and provides a center of rotation for the OAD drive shaft. The *ViperWire* guide wire torquer is an accessory, packaged with the *ViperWire* guide wire, and provides a gripping surface for manipulating the *ViperWire* guide wire.

**Figure 3.** *ViperWire* Guide Wire

#### A. Distal spring tip

#### 1.6 ViperWire Advance and ViperWire Advance with Flex Tip Coronary Guide Wire package contents

- The *ViperWire* guide wire and *ViperWire* guide wire torquer are supplied sterile and are for single-use only. Each package contains:
- Five (5) *ViperWire Advance* guide wires Model GWC-12325LG-FLP or GWC-12325LG-FT (sterile)
  - Five (5) *ViperWire* guide wire torquers (sterile)

#### 1.7 ViperSlide Lubricant

*ViperSlide* Lubricant reduces friction between the OAD drive shaft and the *ViperWire* guide wire. **Note:** Please refer to the *ViperSlide* Lubricant IFU prior to starting the atherectomy procedure.

## 2. Indications for use

The *Diamondback 360* Coronary Orbital Atherectomy System (OAS) is a percutaneous orbital atherectomy system indicated to facilitate stent delivery in patients with coronary artery disease (CAD) who are acceptable candidates for PTCA or stenting due to *de novo*, severely calcified coronary artery lesions.

## 3. Contraindications

Use of the OAS is contraindicated in the following situations:

- The *ViperWire* guide wire cannot pass across the coronary lesion.
- The target lesion is within a bypass graft or stent.
- The patient is not an appropriate candidate for bypass surgery, angioplasty, or atherectomy therapy.
- The patient has angiographic evidence of thrombus.
- The patient has only one open vessel.
- The patient has angiographic evidence of significant dissection at the treatment site.
- Women who are pregnant or children.

## 4. Warnings

- Do **not** use the OAS if the physician does not have experience in coronary angioplasty at their institution.
- Do **not** use the OAS if the physician does not have training on using the OAS. Contact a CSI representative for information on training.
- Do **not** use other commercially-available guide wires with the OAD. Only use the Model GWC-12325LG-FLP *ViperWire Advance* Coronary Guide Wire or GWC-12325LG-FT *ViperWire Advance* Coronary Guide Wire with Flex Tip with the coronary OAD. The *ViperWire* guide wire is designed for use with all coronary OAD crown and shaft configurations.
- **Never** operate the OAD without normal saline and *ViperSlide* Lubricant. Continually flowing saline and *ViperSlide* Lubricant is required for cooling and lubricating the OAD during use in order to avoid overheating and permanent damage to the device and possible injury to the patient.

- Do **not** use the OAD or the *ViperWire* guide wire if their sterile package barriers are compromised or damaged.
- Do **not** use device during spasm of the vessel.
- Do **not** re-sterilize or re-use the OAD. If the OAD is re-sterilized or re-used, the OAD may not function properly potentially leading to serious infection and patient harm and/or death.
- Do **not** re-sterilize or re-use the *ViperWire* guide wire or the guide wire torquer. If the *ViperWire* guide wire or torquer is re-sterilized or re-used, the guide wire may not function properly potentially leading to serious infection and patient harm and/or death.
- **Never** use force to advance the spinning crown as vessel perforation may occur. If any resistance to crown travel is felt, reposition the crown away from the lesion, immediately stop treatment, and use fluoroscopy to assess the vessel for any complications. If it is confirmed there are no complications, reposition the device and advance and retract at a targeted rate of 1 to 3 mm per second.
- Use fluoroscopy to monitor and maintain spacing between the driveshaft and guide wire spring tip throughout the procedure. Always keep more than 5mm of spacing between the distal end of the OAD driveshaft and the proximal end of the guide wire spring tip. If the distance between the driveshaft tip and the *ViperWire* guide wire spring tip is insufficient, the driveshaft tip may contact the guide wire spring tip and result in dislodging the guide wire spring tip.
- Immediately stop using any OAS component should mechanical failure of any component occur before or during the atherectomy procedure. Using damaged components may result in OAS malfunction or patient injury.
- Immediately stop use of the OAD if the device stalls. Review for complications if a stall condition occurs. Do not change to high speed if device stalls on low speed.

**Note:** If a stall occurs, the On/Off button is inactive for five seconds. If the On/Off button is pressed during this five second lockout period, the lockout period will begin again.

- Initial treatment for each lesion must start at low speed.
- Do **not** continue treatment if the wire or the device becomes subintimal.
- Do **not** operate the OAD if there is a bend, kink, or tight loop in the *ViperWire* guide wire. A bend, kink, or tight loop in the *ViperWire* guide wire may cause damage to and malfunctioning of the device during use.
- Performing treatment in excessively tortuous or angulated vessels or bifurcations may result in vessel damage or device failure requiring retrieval.
- Always keep the crown advancing or retracting, while it is orbiting, by continually moving the crown advancer knob to ensure corresponding (1:1) movement between the crown advancer knob and the orbiting crown.
- Do **not** start or stop orbiting of the crown when tight in a lesion.
- Once the OAD has reached full speed (as indicated by a stable pitch) continue to maintain a targeted travel rate of 1 to 3 mm per second, and do not exceed 10 mm per second. If the orbiting crown remains in one location it may lead to vessel damage.
- **Maximum total treatment time should not exceed 5 minutes.** If maximum total treatment time exceeds 5 minutes, the OAD shaft, crown, and *ViperWire* guide wire may begin to exhibit signs of wear and result in a device malfunction and possible injury to patient. A team member should track run time during use to verify total run time is not exceeded.
- Do **not** advance or retract the orbiting crown by advancing the OAD sheath or handle. Buckling of the *ViperWire* guide wire may occur resulting in vessel perforation or vascular trauma. Always advance the orbiting crown by using the crown advancer knob.
- Do **not** inject contrast solution into the OAD injection port. Device failure or patient harm may occur.
- Do **not** attempt aspiration through the OAD or saline line while placed within the body. If saline is pulled out through the OAD or saline line, air may enter the system.
- If air is noticed in the system while the OAD is within the body, discontinue treatment by pressing the OAS Pump power button and carefully remove the OAD driveshaft and crown from the introducer sheath/guide catheter.
- Do **not** allow body parts or clothing to come into contact with spinning components as the OAD orbits at very high speeds. Physical injury to the user or entanglement of clothing with the crown may occur.
- The OAS was only evaluated in severely calcified lesions; therefore the scientific evidence to support use of the OAS to treat other types of lesions/patients is limited.
- Do **not** spin the crown in *SlideAssist*, with the guide wire brake lever in the unlocked position, without first securing the guide wire by holding

it with fingers or by using the guide wire torquer. If using the guide wire torquer, ensure that it is securely fastened to the guide wire before starting to spin the crown. Failure to secure the guide wire when the brake is unlocked could allow the guide wire to spin while in *SlideAssist* mode which may result in patient harm.

## 5. Precautions

- Do **not** use the OAD if there is damage to the OAD package or if the OAD has reached its shelf life expiration date.
- If using an adjustable hemostasis valve with the guide catheter, close the hemostasis valve to minimize blood loss from around the guide catheter while still allowing the OAD sheath to slide through the hemostasis valve. Avoid excessive tightening of the hemostasis valve to prevent damaging the OAD catheter sheath. When inserting or removing the OAD crown or drive shaft through the hemostasis valve, use care not to deform the drive shaft.
- If crown and crown advancer knob movements are not moving correspondingly with one another (1:1 movement), retract and re-advance the crown into the lesion using a travel rate between 1 to 3 mm per second. Repeat retracting and advancing the crown into the lesion until crown to crown advancer knob movement correspondence is observed. If the knob and the crown are not moving together, the crown may be driven into the lesion with too much force and may result in the crown springing forward on exiting the lesion.
- Follow standard institution atherectomy policies and procedures, including those related to anticoagulation, channel blockers, and vasodilator therapy.
- Ensure fluoroscopy provides adequate visualization of the OAS system. Always use fluoroscopy to monitor the guide wire spring tip and driveshaft positions at all times throughout the procedure. If wire movement occurs, it is suggested to reposition the guide wire before advancing the device or continuing treatment.
- A temporary pacing lead may be necessary when treating lesions in the right coronary and circumflex arteries due to the possible occurrence of electrophysiological alternations.
- The risk of the occurrence of a dissection or perforation is increased in severely calcified lesions undergoing percutaneous treatment; therefore, on-site surgical back-up should be included as a clinical consideration.
- Do **not** kink or crush the saline tubing as this will reduce the flow of saline and *ViperSlide* Lubricant to the OAD.
- Continually monitor and check the saline tubing and connections for leaks during the procedure.
- Do **not** spin the crown while advancing or retracting the crown within a guide catheter or tuohy. Damage to the guide catheter, tuohy, and/or OAD may occur.
- Ensure the OAD strain relief remains straight during atherectomy treatment. If the OAD strain relief does not remain straight, the shaft/sheath can kink.
- Do **not** sterilize the OAS pump. Sterilizing will damage the OAS pump. The OAS pump is intended to be used and maintained outside of the sterile field. See Section 10.3 for instructions on cleaning and disinfecting the OAS pump.
- Do **not** allow fluid to leak onto electrical connections of the OAS pump.
- Do **not** spin the crown without a seated and supportive guide catheter
- When treating from a larger lumen to a smaller lumen, make sure the guide catheter is coaxial and that the tip of the OAD has entered the coronary artery to control the initial orbit before engaging the crown; engage the OAD tip into the tight stenosis until low speed has reached its treatment potential prior to initiating treatment with high speed.
- To relieve compression in the driveshaft, lock the crown advancer knob at 1cm from the full back position, advance device over wire to a position proximal from the lesion, deploy the guide wire brake, then unlock the crown advancer knob and move it fully proximal. If the OAD is started with existing compression in the driveshaft it may result in the crown springing forward.
- Do **not** flip contents of tray into sterile field as damage may occur. Components within tray must be carefully removed and placed into sterile field to avoid damage.
- Ejection fractions less than 25% have not been studied, use with low ejection fractions may require additional precautions due to compromised heart function.

## 6. OAS component storage and handling

### 6.1 Storage

Store all OAS components at room temperature in a clean environment away from magnets and sources of electromagnetic interference (EMI). Do not store *ViperSlide* lubricant above 25°C (77°F). Do not freeze *ViperSlide* lubricant. Refer to the *ViperSlide* lubricant IFU prior to starting the atherectomy procedure.

## 6.2 Handling

- All system components are intended to be used in typical operating room/catheterization laboratory environments.
- Additional coronary OAS components should be on hand in the event of damage to any of the components or to component packaging.
- Do not reuse or resterilize the OAD, *ViperWire* guide wire, guide wire torquer, or *ViperSlide* Lubricant as these components are designed for single-use only.
- Do not use the OAD or the *ViperWire* guide wire if their sterile package barriers are compromised or damaged.
- Do not use the OAD or OAS pump if either of them were dropped onto a hard surface, from a height at or greater than 12 in (30 cm), as the OAD or OAS pump may be damaged and may fail to operate properly.
- Do not use any OAS components after their use-by date.
- Do not use *ViperSlide* lubricant if it is exposed to temperatures outside the range indicated on the package labels.
- **Caution:** Do not allow fluid to leak onto electrical connections of the OAS pump.

## 7. Adverse events

Potential adverse events that may occur and/or require intervention include, but are not limited to:

- Allergic reaction to medication/media/device components
- Aneurysm
- Angina (ischemic chest pain)
- Arrhythmias
- Arteriovenous fistula
- Bleeding
- Bruising/hematoma
- Cardiac/cardiopulmonary arrest
- Cardiac/pericardial tamponade
- Cerebrovascular accident (CVA)
- Death
- Embolization, distal (air, tissue, thrombus, device)
- Emergent coronary artery bypass graft surgery (CABG)
- Failure to deliver the system to the intended locations
- Fever
- Heart failure/dysfunction
- Hemorrhage, requiring transfusion
- Hypotension/hypertension
- Infection
- Myocardial infarction
- Pain
- Pericardial effusion
- Pseudoaneurysm
- Restenosis of treated segment leading to revascularization
- Renal insufficiency/failure
- Shock (cardiogenic, hypovolemic)
- Slow flow or no reflow phenomenon
- Stroke
- Thrombus
- Vessel closure, abrupt
- Vessel injury, requiring surgical repair
- Vessel dissection, perforation, rupture, or spasm
- Vessel occlusion

## 8. Clinical trial summary

Clinical trial summary information and data are available at [www.csi360.com](http://www.csi360.com).

## 9. Equipment, set up, and test

### 9.1. Equipment

- In addition to OAS components, equip the operating room with the following:
- Guide catheter - see Table 1 for guide catheter size recommendations
  - Standard IV pole with five wheels and a 20 inch diameter base
  - 1000 mL bag of sterile normal saline
  - Fluoroscopic imaging equipment
  - Standard 100-240V (50/60Hz) hospital grade, electrical wall outlet
  - Other equipment, as needed, for interventional procedures

### 9.2. OAS Pump set up

1. Use the IV pole screw clamp to attach the OAS pump to a standard IV pole making sure to attach the OAS pump to the IV pole at a distance not greater than 60 in (153.0 cm) from the floor to the top edge of the OAS pump.
2. Hang the low saline level sensor and cord, by the closed loop, from the horizontal arm of the standard IV pole.
3. Plug the low saline level sensor connector into the back of the OAS pump (Figure 4).

**Figure 4.** Plug in the low saline level sensor

4. Verify that the power cord is connected to the back of the OAS pump.
5. Insert the other end of the power cord into the electrical wall outlet.

### 9.3. Preparing the ViperWire Advance or ViperWire Advance with Flex Tip Coronary Guide Wire

1. While using sterile techniques, open the ViperWire guide wire packaging pouch and remove the packaging tube.
2. Remove the ViperWire guide wire and ViperWire guide wire torquer from the packaging tube as follows:

- a. Locate the proximal ViperWire guide wire retainer on the inside of the packaging tube. This will expose the proximal end of the ViperWire guide wire.
- b. Removing the ViperWire guide wire – distal end first: advance the proximal end of the ViperWire guide wire into the packaging tube. This will expose the distal end and spring tip of the ViperWire guide wire. Grasp the exposed distal end of the ViperWire guide wire and gently pull the ViperWire guide wire out of the packaging tube. Use care to not stretch or damage the spring distal tip while removing the ViperWire guide wire from the packaging tube.
- c. Removing the ViperWire guide wire – proximal end first: Grasp the exposed proximal end of the ViperWire guide wire and gently pull the ViperWire guide wire out of the packaging tube. Use care to not stretch or damage the spring distal tip while removing the ViperWire guide wire from the packaging tube.
- d. If use of the guide wire torquer is desired, remove the torquer from its packaging. Attach the torquer to the guide wire by holding the distal end of the torquer and rotating the proximal end counterclockwise to tighten. The guide wire is manipulated by firmly gripping the attached torquer. The torquer can be repositioned as necessary.

### 9.4. Initializing the atherectomy procedure

1. Gain vessel access using the physician's preferred methodology.
2. Access the treatment site with a guide catheter.
3. Use angiography to locate, visualize, and evaluate the coronary artery lesion.
4. Approach and cross the lesion, with the ViperWire guide wire, using the physician's preferred methodology.

### 9.5. Preparing the bag of saline and ViperSlide Lubricant

1. Ensure that the OAS pump is powered off by pressing the **Master Power** switch on the back of the OAS pump to off and ensure that no LEDs are illuminated on the OAS pump panel (Figure 5).
2. Prepare a full 1000 mL bag of normal saline solution with ViperSlide Lubricant. Refer to the ViperSlide Lubricant Instructions for Use for lubricant preparation instructions.
3. Hang the prepared saline bag with ViperSlide Lubricant from the low saline level sensor on the standard IV pole.

**Caution:** Do not use glass bottles for the saline solution with ViperSlide Lubricant or hang multiple saline bags from the low saline level sensor as this will disable the Low Saline Information signal.


Figure 5. OAS pump control panel

- A. Low saline red LED indicator
- B. Prime button
- C. Start button and green LED indicator
- D. Status yellow LED indicator

### 9.6. Connecting the OAD to the OAS pump

Remove the sterile saline tubing from the OAD package and pass the saline bag spike end of the saline tubing out of the sterile field. Connect the other end of the saline tubing luer to the OAD luer. Additionally, pass the OAD power cord out of the sterile field.

Perform the following:

1. Connect the saline tubing to the saline bag with ViperSlide Lubricant using standard institution procedures.
2. Open the door, located on the front of the OAS pump, by rotating the door in the direction of the arrow (Figure 6).
3. Place the saline tubing in between the saline tubing positioners that are affixed to the saline tubing, into the top and bottom saline tubing V-guides (Figure 6).


Figure 6. Placing the saline tubing within the OAS pump

- A. Saline tubing positioners
  - B. Saline tubing
  - C. V-guides
  - D. Pump door
4. While closing the door, verify that there is no pinching of the saline tubing and ensure that there is slack in the saline tubing between the OAS pump and saline bag with ViperSlide Lubricant.
  5. Verify that the saline tubing is properly inserted into the saline tubing V-guides and that there are no kinks or damage to the saline tubing.
  6. Press the **Master Power** switch, on the back of the OAS pump, and verify that the red or yellow LED is illuminated on the OAS pump control panel.
  7. Connect the OAD power cord to the OAS pump (Figure 7), the OAS pump, and verify that the red or yellow LED is illuminated on the OAS pump control panel.

Figure 7. Connect the OAD power cord to the OAS pump

8. Remove the drive shaft from the dispenser coil.
9. Push the crown advancer knob fully proximal, away from the nose of the handle.
10. Purge air from the OAD and the saline tubing as follows:
  - a. Verify that the saline tubing is connected to the OAD.
  - b. Press the green **Start** button on the OAS pump control panel to start saline flowing through the saline tubing. Verify that the green LED illuminates.
  - c. Press and hold the **Prime** button on the OAS pump control panel to purge air from the saline tubing. Continually pressing the **Prime** button will pump saline through the tubing at an increasing flow rate. Releasing the **Prime** button will decrease the flow to the low flow rate after two seconds.
  - d. Verify that saline is exiting from the OAD sheath near the crown.
  - e. Continue priming to ensure there are no air bubbles within the saline tubing and use standard hospital procedures to aspirate or purge air from the lines.

**Warning:** Do not attempt aspiration through the OAD or saline line while placed within the body.

- f. After verifying there are no air bubbles within the saline tubing, discontinue priming.

### 9.7. Testing the OAD

#### 9.7.1. Testing OAD crown advancement

Before inserting any portion of the OAD into the body, ensure that axial movement of the OAD crown advancer knob will produce smooth travel of the crown.

**Caution:** Do not spin the crown during this test.

1. Ensure that the crown advancer knob is in the unlocked position as this will allow free axial travel of the crown advancer knob.
2. While visually monitoring the crown, slowly move the crown advancer knob in a back and forth motion. The maximum travel of the crown advancer knob, and the corresponding maximum travel of the shaft tip, is 7.5 cm.

#### 9.7.2. Testing the OAD crown rotation

1. Verify the crown advancer knob is fully proximal, away from the nose of the handle, and release the guide wire brake before threading the ViperWire guide wire through the OAD drive shaft.
2. Grasp the proximal end of the ViperWire guide wire and thread the ViperWire guide wire through the opening in the OAD drive shaft distal tip.

**Warning:** Do not operate the OAD if there is a bend, kink, or tight loop in the ViperWire guide wire. A bend, kink, or tight loop in the ViperWire guide wire may cause damage to and malfunctioning of the OAD during use.

3. Continue feeding the ViperWire guide wire into the OAD drive shaft until the guide wire appears at the rear of the OAD.

**Caution:** If using an adjustable hemostasis valve with the guide catheter, close the hemostasis valve to minimize blood loss from around the guide catheter while still allowing the OAD sheath to slide through the hemostasis valve. Avoid excessive tightening of the hemostasis valve to prevent damaging the OAD catheter sheath. When inserting or removing the OAD crown or drive shaft through the hemostasis valve, use care not to deform the drive shaft.

4. Lock the ViperWire guide wire in place by pressing down on the guide wire brake lever as the crown will not spin if the guide wire brake is unlocked.

5. Verify that saline is still flowing freely out of the saline sheath tip. Verify that the saline tubing is properly connected to the saline bag, that the saline tubing routes correctly through the saline tubing guides, and that the saline tubing is properly connected to the OAD.
6. Press and release the **On/Off** button located on top of the crown advancer knob to activate crown rotation. The OAD is preset to low speed, and the illuminated LED on the OAD will indicate that the OAD is operating at low speed.
7. Check that the flow of saline is increasing and that the shaft and crown are beginning to spin.
8. Immediately press and release the **On/Off** button to stop the shaft and crown from spinning and to complete the test.

### 10. OAS Directions Use

#### 10.1 Performing the atherectomy procedure

1. Ensure that the OAD guide wire brake lever is open (in the up position).
2. Lock the crown advancer knob at 1 cm from the fully proximal position by rotating the crown advancer knob 90 degrees in either direction.
3. Advance the OAD drive shaft over the ViperWire guide wire and through the hemostasis valve while keeping ViperWire guide wire placement stationary.
4. While using fluoroscopy, gently advance the OAD crown over the ViperWire guide wire to a position outside the lesion in a non-occlusive location. Verify that the OAD distal tip is not within the lesion when the crown and drive shaft begin to spin.

**(Optional)** Use the *GlideAssist* feature to facilitate advancing the OAD crown over the ViperWire guide wire.

**Warning:** Spinning the crown using *GlideAssist* can be done with the OAD guide wire brake lever in either the locked or unlocked position. If using *GlideAssist* with the guide wire brake in the unlocked position, the guide wire must be held using either fingers or the guide wire torquer. If using the guide wire torquer, ensure that it is securely fastened to the guide wire before starting to spin using *GlideAssist*.

- a. Enable *GlideAssist* by pressing and holding the low speed button. Release the button once the low speed light begins to slowly blink. The slowly blinking light indicates *GlideAssist* is enabled.
- b. Ensure that the guide wire is secure by locking the guide wire brake or by holding the guide wire with either fingers or the guide wire torquer.
- c. Spin by pressing and releasing the **On/Off** button. The crown will begin spinning and the low speed light will rapidly blink indicating that the crown is spinning in *GlideAssist* mode.
- d. Stop spinning by pressing and releasing the **On/Off** button. The low speed light will slowly blink indicating that the OAD is no longer spinning but continues to be in *GlideAssist* mode.
- e. Disable *GlideAssist* by pressing and immediately releasing either the low or the high speed button while the crown is not spinning. The low speed light will stop blinking, yet remains illuminated indicating that the OAD is now in treatment mode.

**Note:** If the brake configuration is changed from either the locked or unlocked position while spinning in *GlideAssist* mode, the crown will automatically stop spinning yet the OAD will remain in *GlideAssist* mode.

5. Inject contrast medium, through a port in the hemostasis valve, to verify that the size of the crown is compatible with the treatment area diameter (see Appendix B). Verify that contrast media injections are not above 400 PSI.
6. Verify that the ViperWire guide wire spring tip is distal to the lesion and is not in danger of coming in contact with the spinning crown and drive shaft tip.

**Caution:** Maintain at least 5 mm between the proximal end of the ViperWire guide wire spring tip and the OAD drive shaft tip to prevent contact of the drive shaft tip with the guide wire spring tip. Further advance the ViperWire guide wire, as necessary, to maintain the 5 mm minimum distance.

7. Push down on the ViperWire guide wire brake lever to engage the brake. The crown will not spin if the guide wire brake is not locked.
8. Unlock and push the crown advancer knob to the fully proximal position to relieve any compression in the driveshaft.
9. Press and release the **On/Off** button on top of the crown advancer knob to begin spinning the crown. The OAD is preset to low speed, and the illuminated LED on the OAD will indicate that the OAD is operating at low speed.

**Warning:** Initial treatment for each lesion must start at low speed.

**Caution:** Continually monitor the saline fluid levels during the procedure. Continual infusion of saline and ViperSlide Lubricant is critical for safe coronary OAS operation.

10. Audibly verify that the OAD drive shaft and crown are orbiting at a stable speed as indicated by the OAD frequency (pitch) stabilizing following the 1.5 second ramp up in speed.
11. Slowly advance and retract the crown advancer knob to begin atherectomy of the lesion at a targeted travel rate of 1 to 3 mm per second, and do not exceed 10 mm per second. During initial passes of a lesion, or when treating a tight lesion, advance and retract slowly until crown to crown advancer knob movement correspondence if observed (1:1 movement). The device works both antegrade and retrograde, avoid rapid movements while the device is orbiting. Using fluoroscopy, continually verify that the crown and the crown advancer knob are moving correspondingly with one another. Ensure that the OAD remains horizontal during the procedure to minimize saline leakage from the OAD handle.

12. Using a series of intermittent treatment intervals and rest periods, slide the crown advancer knob to move the crown back and forth across the lesion always returning to a non-occlusive position within the vessel when the interval set is complete.

**Warning:** Once OAD has reached full speed (as indicated by a stable pitch) do not allow the orbiting crown to remain in one location as it may lead to vessel damage. Continue to maintain a targeted travel rate of 1 to 3 mm per second, and do not exceed 10 mm per second.

For every 30 seconds of treatment, a rest period of equal time is recommended with a maximum treatment time of 5 minutes per OAD. A team member should track individual treatment times to verify treatment does not exceed 30 seconds. A team member should also track that corresponding rest periods are following individual treatment times as well. The OAS pump will emit a beep after every 25 second interval of treatment time.

**Warning: Maximum total treatment time should not exceed 5 minutes.** If maximum total treatment time exceeds 5 minutes, the OAD shaft, crown, and ViperWire guide wire may begin to exhibit signs of wear and result in a device malfunction and possible injury to patient. A team member should track run time during use to verify total run time is not exceeded.

#### 10.1.1. Replacing the bag of saline and ViperSlide Lubricant

The low saline level sensor triggers an audible information signal every 5 seconds for a total of 30 seconds if there is less than 200 mL (+/- 100 mL) remaining in the bag of saline and ViperSlide lubricant during a treatment period. If the low saline level sensor triggers during a rest period, only the red low saline LED is illuminated. Perform the following to replace the bag of saline and ViperSlide Lubricant:

1. Ensure that the pump is stopped by pressing the green **Start** button on the OAS pump control panel and verify that the green LED, on the OAS pump control panel, is not illuminated.
2. Prepare a new 1000 mL bag of normal saline solution with ViperSlide Lubricant. Refer to the ViperSlide Lubricant Instructions for Use for lubricant preparation instructions.
3. Remove the empty bag of saline and ViperSlide Lubricant from the low saline level sensor on the IV pole.
4. Hang the new bag of saline and ViperSlide Lubricant from the low saline level sensor on the standard IV pole.

**Caution:** Do not use glass bottles for the saline solution with ViperSlide Lubricant or hang multiple saline bags from the low saline level sensor as this will disable the Low Saline Level Sensor.

5. Remove the bag spike from the empty bag of saline and ViperSlide Lubricant and spike the new bag of saline and ViperSlide Lubricant.
6. Power on the OAS pump by pressing the green **Start** button on the OAS pump control panel.
7. Ensure that no air was introduced into the saline tubing.

### 10.1.2. Replacing the OAD

If the OAD needs replacing, perform the following:

1. Stop the spinning crown and drive shaft by pressing and releasing the **On/Off** button on top of the crown advancer knob.
2. Disconnect the OAD power cord from the OAS pump.
3. Leave the guide catheter and the *ViperWire* guide wire in place, release the guide wire brake on the OAD, and retract the OAD sheath and drive shaft, from the guide catheter, while monitoring and maintaining the current guide wire position.
4. Power off the OAS pump by pressing the green **Start** button on the OAS control panel to stop saline from flowing through the saline tubing and verify that the green LED, on the OAS pump control panel, is not illuminated.
5. Disconnect the saline tubing from the OAD currently in use and set aside for use with the replacement OAD.
6. Obtain a new replacement OAD and remove the new replacement OAD from the package.
7. First, attach the existing saline tubing to the new replacement OAD, then connect the new replacement OAD power cord to the OAS pump.
8. Press the green **Start** button on the OAS pump control panel to start the saline flowing through the saline tubing and verify that the green LED illuminates.
9. Purge the air from the OAD. Refer to Section 9.6, step 10.
10. Load the new replacement OAD drive shaft over the existing *ViperWire* guide wire.
11. Test the OAD crown advancement per the instructions in Section 9.7.1.
12. Test the OAD crown rotation per the instructions in Section 9.7.2.

### 10.2. Completing the atherectomy procedure

To complete the atherectomy procedure, perform the following:

1. While the crown is spinning, retract the crown and drive shaft proximal to the lesion at a targeted travel rate of 1 to 3 mm per second, and do not exceed 10 mm per second.
2. Stop the OAD crown and drive shaft rotation by pressing and releasing the **On/Off** button on top of the crown advancer knob.
3. Release the guide wire brake on the OAD and carefully remove the OAD drive shaft and crown from the guide catheter.

**(Optional):** Use the *SlideAssist* feature to facilitate retracting the OAD crown over the guide wire.

**Caution:** Do not spin while retracting the crown within a guide catheter or tuohy. Damage to the guide catheter, tuohy, and/or OAD may occur.

4. Press the green **Start** button on the OAS pump control panel to stop saline from flowing through the saline tubing and verify that the green LED is not illuminated. Turn off the OAS pump at the **Master Power** switch on the back of the OAS pump.
5. Remove and dispose of the *ViperWire* guide wire and guide catheter according to standard hospital procedures.
6. Treat the puncture site according to standard interventional procedure protocol.

**Warning:** The OAD, *ViperWire* guide wire and *ViperSlide* Lubricant are designed for single patient use only and should not be reused or re-sterilized. The saline tubing and partially-used bag of saline and *ViperSlide* Lubricant are designed for single patient use only and should not be stored or reused. Discard devices according to hospital guidelines at the end of the procedure.

**Note: Disposal of the OAD:** The OAD is designed for single use and should not be reused or re-sterilized. Discard the OAD and saline tubing according to standard hospital protocol.

Certain states may have additional requirements for disposal of certain batteries. Please verify your state's requirements for disposal of CR lithium coin cells (CR2032 or equivalent) prior to disposal.

**California, USA Only:** Perchlorate Material – special handling may apply. See [www.dtsc.ca.gov/hazardouswaste/perchlorate](http://www.dtsc.ca.gov/hazardouswaste/perchlorate)

### 10.3. Maintaining the OAS Pump

The OAS pump does not require routine maintenance, periodic maintenance, or calibration. *CSI* recommends inspection of the OAS pump in accordance with the hospital's standard biomedical engineering department protocol. The pump has been designed to function for 875 hours minimum, with 350 hours of minimum OAD use, which equates to 5 years. Contact *CSI* Customer Service if there are questions about the OAS pump function or performance.

### 10.3.1. Cleaning the OAS Pump

Clean the OAS pump immediately after each use by following the steps below:

**Caution:** Ensure that the OAS pump is powered off at the **Master Power** switch on the back of the OAS pump and disconnect the OAS pump from wall power before cleaning the OAS pump.

**Caution:** Do not immerse the OAS pump into fluids. Do not use solvents or abrasive cleaners to clean the OAS pump as these may damage the OAS pump and OAS pump components.

**Caution:** Completely dry the OAS pump before reconnecting the OAS pump to wall power and powering on the OAS pump.

1. Prepare an enzymatic detergent, such as Enzol, per manufacturer's directions.
2. Thoroughly wipe down the OAS pump, using a clean soft cloth that has been dampened with the prepared detergent, until all visible soil is removed.
3. Thoroughly rinse the OAS pump, using a clean soft cloth that has been dampened with lukewarm tap water.
4. Dry the OAS pump using a clean, soft cloth; and, if available, filtered, pressurized air at ≤40 psi.

### 10.3.2. Disinfecting the OAS Pump

1. Verify no debris is present after pump has been cleaned and rinsed with enzymatic detergent. Repeat the above cleaning procedure if any debris continues to be visible.
2. Put on a pair of disposable protective gloves. Check the expiration date on container and remove a fresh moist Super Sani-Cloth. Wring excess solution from the wipe, ensuring it is saturated, but not dripping. Discard wipes as they become dry.
3. Disinfect the front pump face, ensuring that all the indicated surfaces are completely covered with solution during the wiping process for a total of 2 ½ minutes to ensure an adequate "dwell" time. Dwell time means the number of minutes that a product must be in contact with the surface, and remain wet, in order to assure proper efficacy, or effectiveness to kill organisms. **Surfaces must not become dry at any point during disinfection.** Refer to the following steps for surfaces to be wiped:
  - a. Open the pump head cover (see Figure 8).
  - b. Thoroughly wipe the edge of the pump head cover all along the closure seam on both sides (see Figures 9 and 10 below).

Figure 8. Open pump head cover


Figure 9. Wipe closure seam

Figure 10. Wipe closure seam

4. Thoroughly wipe the seams and crevices of the pump head cover, around the edges of the lettering, around the button area, and around the power outlet (see Figure 11). Thoroughly wipe all surfaces of the front face of the pump (see Figure 12). Continue to keep all surfaces wet for a minimum of two and a half minutes. Discard the wipes.

Figure 11. Thoroughly wipe these areas



Figure 12. Wipe all surfaces of front face

### 10.3.3. Disinfecting the OAS Pump – users outside the United States

Disinfect the OAS pump after each use by following the steps below:

**Caution:** Ensure that the OAS pump is powered off at the **Master Power** switch on the back of the OAS pump and disconnect the OAS pump from wall power before disinfecting the OAS pump.

**Caution:** Do not immerse the OAS pump into fluids. Do not use solvents or abrasive disinfectants to disinfect the OAS pump as these may damage the OAS pump and OAS pump components.

**Caution:** Completely dry the OAS pump before reconnecting the OAS pump to wall power and powering on the OAS pump.

1. Open a fresh, sterile wipe that is pre-saturated with 70% Isopropyl Alcohol (IPA) or prepare a sterile gauze/wipe by pouring or soaking it with 70% IPA. Wring any excess IPA from the gauze/wipe, ensuring that the gauze/wipe remains saturated, but not dripping, with IPA.
2. Thoroughly wipe all surfaces on the front face of the OAS pump. Concentrate wiping the seams and crevices of the OAS pump head door, around the edges of the lettering, and around the OAS pump control panel. Continue to wipe these surfaces for a minimum of one (1) minute. Discard the gauze/wipe. Repeat step #1.
3. Open the OAS pump head door. Using the prepared gauze/wipe, thoroughly wipe the edge of the OAS pump head door on both sides of the closure seam. Continue to wipe these surfaces for a minimum of one (1) minute. Discard the gauze/wipe. Repeat step #1.
4. Thoroughly wipe all surfaces on the front face of the OAS pump and the edge of the OAS pump head door on both sides of the closure seam. Continue to wipe these surfaces for a minimum of one (1) minute. Discard the gauze/wipe. Repeat step #1.
5. Continue to repeat this wiping process, as many times as necessary, to ensure that all surfaces remain wet with IPA for a minimum of ten (10) minutes.
6. Close the OAS pump door when disinfecting is complete.

### 10.4. Returning OAS components

Contact *CSI* Customer Service for returning OAS components. See the back of this instruction for use for *CSI* contact information.

## 11. Specifications

### 11.1. OAD specifications

Parameter	Value
Nose length	5 mm
Maximum OAD sheath outer diameter	1.45 mm
Minimum guide catheter diameter	6 French with an internal diameter of at least 0.066 in (1.68 mm)
Electrical cable length: OAD to OAS pump	3.4 m (11 ft)
Electrical connector type (device power)	Type CF applied Part – DC barrel (48 V DC)
Fluid connector type	Luer fitting
Tubing length (from saline bag to OAD)	3.7 m (12 ft)

Parameter	Value
Visual alerts	Speed indicators
Sterilization	Ethylene oxide (EtO) cycle
Rotational speeds	5000 rpm 80000 rpm 120000 rpm
Operating life	5 minutes of total therapy time
Water Ingress Protection	IPX1: Protection against water ingress
Approximate saline flow rate	
•Prime button pressed, OAD not spinning	16-36 ml/min
•OAD spinning on low speed	11-34 ml/min
•OAD spinning on high speed	7-29 ml/min
•OAD not spinning, prime button not pressed	12-19 ml/min

## 11.2 OAS Pump Specifications

Parameter	Value
Depth	<30.6 cm (12.0 in)
Height	20.3 cm (8.0 in)
Width	25.4 cm (10.0 in)
Weight	<5.0 kg (11 lbs)
Electrical cable length: OAS pump to electrical outlet	6.1 m (20ft)
Master Fuse	250 V 4A SLOW BLOW 1/4" x 1 - 1/4"
External housing	ABS Plastic
Electrical connector type (Main Power)	Mains Power Plug (100-240V AC @ 50-60 Hz)
Audible information signals	Audible information signal for approximately every 25 sec of OAD spin time.* Audible information signal every 5 sec for a total of 30 sec when the saline level falls below 200 mL during a treatment period.
Visual alerts	Start button Low Saline Information Signal when ≤200 mL (±100 mL) of 1000-mL bag of saline remaining
Storage conditions	Room temperature in a clean environment
Operating conditions	Typical operating room/catheterization laboratory environment (10-30°C)
Operating life	875 hours minimum, with 350 hours of therapy minimum or 5 years
Water Ingress Protection	IPX1: Protection against water ingress

\*Timer resets when spinning of the OAD stops.

## 12. OAS Pump Declaration of Conformity

*CSI* declares that the coronary OAS is in conformity with the requirements of IEC 60601-1. The OAS pump is compatible for use in a standard catheter laboratory environment.

## 13. EMC Declaration

Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Install and use medical electrical equipment according to the EMC information below:

- Do not have portable and/or mobile radio-frequency (RF) communications equipment within close proximity of medical electrical equipment as portable and mobile RF communications equipment can affect medical electrical equipment.
- Ensure that power frequency magnetic fields are at levels characteristic of a typical commercial or hospital environment.
- Under an EMC phenomena the OAS may stop operation, and may require user intervention to recycle the power to resume operation.

The Orbital Atherectomy System has been tested to IEC 60601-1-2. The Orbital Atherectomy System has been tested to Immunity and Emission Test Levels of a Professional Healthcare Facility Environment. The OAS is Group 1 (Therapy ME Equipment and Systems) and therefore must meet CISPR 11 Class A.

**Note:** The Emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment

**Essential Performance of the Coronary Classic OAD Product:**

- System starts, runs, stops at the operator's discretion and recovers in a controlled manner from external upsets.

**Emissions**

Emissions Standard	Test	Compliance Level	Electromagnetic environment - Guidance
IEC 60601-1-2	RF Emissions CISPR 11	Group 1	The orbital atherectomy device does not utilize RF energy for its function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
IEC 60601-1-2	RF Emissions CISPR 11	Class A	The orbital atherectomy device is suitable for use in all locations other than those located in residential environments and those directly connected to a low-voltage power supply network that supplies buildings used for domestic purposes.
IEC 60601-1-2	Harmonics, IEC 61000-3-2	Class A	
IEC 60601-1-2	Flicker, IEC 61000-3-3	Complies	

**Immunity**

Immunity Standard	Test	Test Level	Compliance Level
IEC 60601-1-2	Electrostatic Discharge (ESD), IEC 61000-4-2	±8 kV contact ±2, ±4, ±8, ±15 kV air	±8 kV contact ±2, ±4, ±8, ±15 kV air
IEC 60601-1-2	Radiated RF, IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz
IEC 60601-1-2	Radiated RF, Proximity Fields, IEC 61000-4-3	Tested to levels specified in Table 9 of IEC 60601-1-2:2014	Complies to levels of Table 9
IEC 60601-1-2	Electrical Fast Transient / Burst, IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines
IEC 60601-1-2	Surge, IEC 61000 4-5	± 2 kV line to ground ± 1 kV line to line	± 2 kV line to ground ± 1 kV line to line
IEC 60601-1-2	Conducted Disturbances RF, IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM bands	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM band
IEC 60601-1-2	Power Frequency 50/60 Hz Magnetic Field, IEC 61000-4-8	30 A/m, 50 / 60 Hz	30 A/m, 50 and 60 Hz
IEC 60601-1-2	Voltage Dips and Interruptions, IEC 61000-4-11	100 % dip for 0.5 cycle; 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 100% dip for 1 cycle 30% dip for 25 cycles 100% dip for 5 seconds	100% dip for 0.5 cycle; 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 100% dip for 1 cycle 30% dip for 25 cycles 100% dip for 5 seconds
AIM 7351731	ISO 14223	65 A/m, 134.2 kHz	65 A/m
AIM 7351731	ISO/IEC 14443-3 (Type A)	7.5 A/m, 13.56 MHz	7.5 A/m
AIM 7351731	ISO/IEC 14443-4 (Type B)	7.5 A/m, 13.56 MHz	7.5 A/m
AIM 7351731	ISO/IEC 15693 (ISO 18000-3 Mode 1)	5 A/m, 13.56 MHz	5 A/m
AIM 7351731	ISO 18000-3 Mode 3	12 A/m, 13.56 MHz	12 A/m
AIM 7351731	ISO/IEC 18000-7	3 V/m, 433 MHz	3 V/m
AIM 7351731	ISO/IEC 18000-63 Type C	54 V/m, 860 – 960 MHz	54 V/m
AIM 7351731	ISO/IEC 18000-4 Mode 1	54 V/m, 2.45 GHz	55 V/m

Appendices

**A. OAS troubleshooting**

If issues with the OAD cannot be resolved in each of the situations below, replace the OAD and continue with the procedure. Contact CS/ Customer Service for returning OAS components. See the back of the instructions for use for CS/ contact information.

Issue Number	Issue	Solution
1	The crown stops spinning during the procedure	<ol style="list-style-type: none"> <li>1. Immediately discontinue treatment. Stop the OAD from spinning, but leave the OAS pump running.</li> <li>2. Verify that saline is flowing.</li> <li>3. Check that the OAS pump green LED light is on and that the OAD green LED light is on. If the OAS pump green LED is flashing, press the pump start button twice.</li> <li>4. Check to ensure that the OAS pump power cord is connected to the back of the OAS pump and that the OAD power cord is connected to the OAS pump.</li> <li>5. Check that the OAD guide wire brake lever is in the down/locked position.</li> <li>6. Position the crown outside the lesion in a non-occlusive location. Use fluoroscopy to analyze the situation prior to attempting a low speed spin of the crown.</li> </ol>
2	Blood is backing up into the OAD	<ol style="list-style-type: none"> <li>1. Immediately discontinue treatment. Stop the OAD from spinning, but leave the OAS pump running.</li> <li>2. Verify that the saline tubing is properly connected to the saline bag, that the saline tubing is routed correctly through the OAS pump saline tubing guides, and that the saline tubing is properly connected to the OAD.</li> <li>3. If the saline tubing is properly connected and blood continues to back into the OAD sheath, replace the OAD.</li> </ol>
3	Crown rotational speeds are variable and will not stabilize	<ol style="list-style-type: none"> <li>1. Immediately discontinue treatment.</li> <li>2. Stop the OAD from spinning, but leave the OAS pump running.</li> <li>3. Verify that saline is flowing.</li> <li>4. Verify <i>ViperSlide</i> Lubricant is present in the saline bag. See the <i>ViperSlide</i> Lubricant Instructions for Use for information. Verify that the saline tubing is properly connected to the saline bag, that the saline tubing is routed correctly through the OAS pump saline tubing guides, and that the saline tubing is properly connected to the OAD.</li> <li>5. Verify that the crown advancer knob moves smoothly.</li> <li>6. Position the crown outside the lesion in a non-occlusive location. Using a targeted travel rate of 1 to 3 mm per second, and not exceeding 10 mm per second, continue treatment on low speed.</li> </ol>
4	All speed indicator LEDs on the OAD handle remain illuminated	Immediately discontinue treatment and replace the OAD.
5	The crown is not moving one-to-one with the crown advancer knob	<p>During start up in the vessel:</p> <ol style="list-style-type: none"> <li>1. Verify the Tuohy valve is not over-tightened.</li> <li>2. Verify the crown advancer knob moves smoothly.</li> <li>3. Retract the crown advancer knob until the crown moves with the knob.</li> </ol> <p>While spinning:</p> <ol style="list-style-type: none"> <li>1. Immediately discontinue treatment. Stop the OAD from spinning, but leave the OAS pump running.</li> <li>2. Verify the Tuohy valve is not over-tightened.</li> <li>3. Verify that the crown advancer knob moves smoothly.</li> <li>4. Retract the crown advancer knob until the crown moves with the knob.</li> <li>5. Verify that contrast media injections are not above 400 psi and are not occurring during spinning of the crown.</li> <li>6. Engage and disengage the lesion using a targeted travel rate of 1 to 3 mm per second, not exceeding 10 mm per second, and verify that the crown and crown advancer knob are moving correspondingly with one another (1:1 movement).</li> </ol>
6	The OAD will not exit <i>SlideAssist</i> mode	Press and immediately release the speed button. Do not continuously hold the speed button to exit <i>SlideAssist</i> mode.
7	The OAS pump will not power on and no LEDs are illuminated on the OAS pump control panel	Ensure that the power cord is properly inserted into the power module on the back of the OAS pump and that the power cord is connected to a functioning wall power outlet. Ensure that the <b>Master Power</b> switch, on the back of the OAS pump, is in the on position.
8	The OAS pump will not pump saline	<ol style="list-style-type: none"> <li>1. Ensure that the OAS pump is properly powered on – see Issue number 6.</li> <li>2. Ensure that the saline bag and saline tubing (i.e. bag spike) are properly connected and a sufficient amount of saline is in the saline bag such that the low saline level sensor is not active and the red LED on the OAS pump control panel is not illuminated.</li> <li>3. Ensure that the OAS pump head is properly aligned with the OAS pump base.</li> <li>4. Ensure that the saline tubing is routed correctly through the OAS pump saline tubing guides and that the OAS pump saline tubing door is closed.</li> <li>5. Ensure that the yellow LED is off and the green LED is illuminated. If the green LED is not illuminated, press the green <b>Start</b> button and verify that the yellow LED is off and that the green LED illuminates.</li> <li>6. Ensure that the yellow LED is off and the green LED is illuminated. If the green LED is flashing while the yellow LED is illuminated, press the green <b>Start</b> button twice and verify that the yellow LED is off and that the green LED illuminates.</li> </ol>

Issue Number	Issue	Solution
9	The green (start) LED is illuminated, but the OAS pump does not pump and the yellow LED is illuminated as well	Ensure that the saline tubing is routed correctly through the OAS pump saline tubing guides and that the OAS pump saline tubing door is closed. Contact CS/ Customer Service at the phone number on the back of this instructions for use.
10	The OAS pump was running, but has stopped pumping and the yellow LED is illuminated	Press the <b>Master Power</b> switch, on the back of the OAS pump, to power off the OAS pump. Wait five (5) seconds and press the <b>Master Power</b> switch to power on the OAS pump.
11	The low saline level sensor (red LED) is illuminated	<p>Note: The OAS pump will stop pumping saline and supplying power to the OAD 30 seconds after the low saline level sensor activates while the OAD is spinning, as indicated by an audible information signal every 5 seconds.</p> <ol style="list-style-type: none"> <li>1. If there is less than 200 mL of saline left in the bag of saline and <i>ViperSlide</i> Lubricant, replace the bag with a new 1000 mL bag of normal saline and <i>ViperSlide</i> Lubricant solution.</li> <li>2. Ensure that the bag of saline and <i>ViperSlide</i> Lubricant is hanging freely from the saline bag open hook and that the low saline level sensor cord is properly inserted into the connector on the sensor and the connector on the back of the OAS pump.</li> <li>3. Verify that the red low saline LED on the OAS pump control panel turns off and either the yellow LED or the green LED illuminates. If the yellow LED illuminates, press the <b>Start</b> button on the OAS pump and verify that the green LED illuminates.</li> </ol>
12	The OAS pump will not power on	<p>Note: The OAS pump is designed to stop pumping saline and cut the supply of power to the OAD if the OAS pump encounters internal errors, such as: improper OAS pump pumping speed, problems with the power supply, inadvertent opening of the saline tubing door, etc. If an internal error occurs, the green LED will turn off and the yellow LED will illuminate.</p> <ol style="list-style-type: none"> <li>1. Ensure that the saline tubing door is closed.</li> <li>2. Attempt to restart the OAS pump by pressing the green <b>Start</b> button on the OAS pump control panel and verify that the green LED illuminates.</li> <li>3. If the OAS pump does not restart after completing the above mentioned steps, press the <b>Master Power</b> switch to power off the OAS pump. Wait a few seconds and press the <b>Master Power</b> switch to power on the OAS pump. Verify that the OAS pump powers on.</li> </ol>
13	All OAD speed indicator LEDs blink simultaneously.	<ol style="list-style-type: none"> <li>1. Discontinue treatment.</li> <li>2. Turn off the power module switch on the back of the OAS pump. After a few seconds, turn the power module switch to on and press the green <b>Start</b> button.</li> <li>3. If the OAD LEDs continue to blink simultaneously, replace the OAD.</li> </ol>
14	The crown cannot spin.	<ol style="list-style-type: none"> <li>1. Discontinue treatment.</li> <li>2. Check that the brake is locked.</li> <li>3. Check the OAS pump and device indicator LEDs.</li> <li>4. Check the power cords.</li> <li>5. Verify that the OAS pump and device are receiving power.</li> </ol>
15	There is no saline flow.	<ol style="list-style-type: none"> <li>1. Discontinue or do not initiate treatment.</li> <li>2. Verify that the saline tubing is properly connected and saline is flowing.</li> <li>3. Ensure that saline tubing is installed properly in the roller OAS pump.</li> <li>4. Verify that the green LED on the front of the OAS pump is on after pressing the <b>Power</b> button.</li> <li>5. If the device is properly connected and no flow is noted, replace the saline tubing, bag, and/or device.</li> <li>6. Verify that the OAS pump door is closed.</li> <li>7. Verify that there are no kinks in all saline tubing and the device saline sheath.</li> </ol>
16	No LEDs lit on front panel of the OAS pump.	<ol style="list-style-type: none"> <li>1. Ensure that the OAS pump is properly powered up (see above).</li> <li>2. If none of the LEDs on the front cover light up after powering the OAS pump, contact CS/ customer service.</li> </ol>
17	The crown will not stop spinning	<ol style="list-style-type: none"> <li>1. Press the power button on the OAS pump.</li> <li>2. Unplug the device from the OAS pump.</li> <li>3. Unplug the OAS pump from the power source.</li> </ol>
18	All three LEDs on the front panel of the OAS pump remain illuminated.	<ol style="list-style-type: none"> <li>1. Press the <b>Master Power</b> switch, on the back of the OAS pump, to power off the OAS pump. Wait a few seconds and press the <b>Master Power</b> switch to power on the pump.</li> <li>2. Contact CS/ Customer Service at the phone number on the back of this instructions for use.</li> </ol>
19	After OAS pump power up, all three LEDs on the front panel of the OAS pump blink three times and there is an audible notification signal three times.	Contact CS/ Customer Service at the phone number on the back of this instructions for use.

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### B. Maximum orbit and result lumen diameter

The following table shows the maximum orbit and resulting lumen diameter for the classic crown size at incremental rotational speeds after approximately 5 minutes of orbital treatment time. The table compares resulting lumens for travel rates of 1 mm per second and 10 mm per second.

**Note:** A pass is defined as once out and back across the lesion. Orbit data presented are based on a 6 cm pass distance.

#### Maximum Lumen Diameters after Approx. 5 Minutes for Different Travel Rates

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### C. Orbit performance

The following chart demonstrates typical orbit diameter vs. duration of operation (as measured in simulated calcified lesions.) This chart is for reference only. Actual orbit performance may vary.

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