

THROMBUS MANAGEMENT

PRODUCT CODE	NAME	OUTER DIAMETER AND CROSSING PROFILE **	INNER DIAMETER**	SHEATH COMPATIBILITY	WORKING LENGTH
A5-055	ASPIRE MAX 5 Aspiration Catheter System* 55cm	5F (0.066")	4.1F (0.054")	5F	55cm
A5-135	ASPIRE MAX 5 Aspiration Catheter System* 135cm	5F (0.066")	4.1F (0.054")	5F	135cm
X5-055	MAX 5 Aspiration Catheter 55cm Only	5F (0.066")	4.1F (0.054")	5F	55cm
X5-135	MAX 5 Aspiration Catheter 135cm Only	5F (0.066")	4.1F (0.054")	5F	135cm
A6-055	ASPIRE MAX 6 Aspiration Catheter System* 55cm	5.9F (0.076")	4.9F (0.064")	6F	55cm
A6-135	ASPIRE MAX 6 Aspiration Catheter System* 135cm	5.9F (0.076")	4.9F (0.064")	6F	135cm
X6-055	MAX 6 Aspiration Catheter 55cm Only	5.9F (0.076")	4.9F (0.064")	6F	55cm
X6-135	MAX 6 Aspiration Catheter 135cm Only	5.9F (0.076")	4.9F (0.064")	6F	135cm
ASP-30	ASPIRE ASPIRATOR 30ml Only				

*System includes ASPIRE Aspirator, MAX Aspiration Catheter, and Connection Tubing. Requires drain bag and Y-adapter. Single use only.

**Reference: 4F = 0.053" 5F = 0.066" 6F = 0.079"