


**Features:**

- Precise tip deflection adjustment for fast access
- Tip position indicator showing neutral or deflected curve
- Superior curve retention for reliable catheter placement
- Enhanced radiopaque tip visibility enabling accurate tip positioning
- Proprietary hydrophobic coated sheath technology, enabling smooth access
- Reliable Flexsteer™ reinforced braid technology provides exceptional control, kink resistance, while maintaining optimal flexibility
- Ergonomic handle with clear imprint of the curve deflection, french size and sheath length
- Reliable SureSeal™ hemostatic valve technology to minimize air embolism and blood loss even under extreme conditions

**ORDERING INFORMATION**



**INTERVENTIONAL USE - ORDERING INFORMATION**

Model # Set <i>(One per box)</i>	A1 Sheath Inner Diameter	A2 Usable Sheath Length	B Usable Dilator Length	A3 Curve Bend
<b>STANDARD</b>				
DST0654509	● 6.5 F	45 cm	63 cm	9 mm
DST0655509	● 6.5 F	55 cm	73 cm	9 mm
DST0659009	● 6.5 F	90 cm	108 cm	9 mm
DST0704509	● 7 F	45 cm	63 cm	9 mm
DST0705509	● 7 F	55 cm	73 cm	9 mm
DST0705517	● 7 F	55 cm	73 cm	17 mm
DST0709009	● 7 F	90 cm	108 cm	9 mm
<b>TD CURVE</b>				
DST0654509TD	● 6.5 F	45 cm	63 cm	9 mm
DST0655509TD	● 6.5 F	55 cm	73 cm	9 mm
DST0704509TD	● 7 F	45 cm	63 cm	9 mm
DST0705509TD	● 7 F	55 cm	73 cm	9 mm

**TRANSSEPTAL USE - ORDERING INFORMATION**

Model # Set <i>(One per box)</i>	A1 Sheath Inner Diameter	A2 Usable Sheath Length	B Usable Dilator Length	A3 Curve Bend	C Compatible Transseptal Needle
<b>STANDARD</b>					
DST0856717	● 8.5 F	67 cm	85 cm	17 mm	89 cm
DST0857317	● 8.5 F	73 cm	91 cm	17 mm	98 cm
DST0856722	● 8.5 F	67 cm	85 cm	22 mm	89 cm
DST0857322	● 8.5 F	73 cm	91 cm	22 mm	98 cm
DST0856750	● 8.5 F	67 cm	85 cm	50 mm	89 cm
DST0857350	● 8.5 F	73 cm	91 cm	50 mm	98 cm

*Guiding sheaths are supplied sterile. For single use only. Available bulk upon request. Standard lengths are noted above. Other lengths available upon request.*



P: 727.937.2511  
F: 727.934.9835  
E: sales@oscor.com  
www.oscor.com

**CAUTION:** Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.





AUG 2014



*Destino™* **TWIST**  
DEFLECTABLE STEERABLE GUIDING SHEATH


The *Destino™ TWIST* is a steerable guiding sheath that offers a deflectable tip with an ergonomic handle. The Destino™ Twist is the ultimate tool for gaining access to the most difficult to reach sites, while maintaining hemostasis with Oscor SureSeal™ technology. The Destino™ Twist is **YOUR ONE SHEATH SOLUTION** intended to facilitate the intracardiac, renal, and peripheral placements of diagnostic and therapeutic devices.




### INTERVENTIONAL USE




#### SHAFT MODELS



#### SHAFT DIAMETERS



#### TIP DEFLECTION




### TRANSEPTAL USE



#### SHAFT MODELS



#### SHAFT DIAMETERS



#### TIP DEFLECTION

