



REACH

**NEW
PRODUCT**



DESTINO™

REACH

BI-DIRECTIONAL DEFLECTABLE GUIDING SHEATH





FEATURES:




- Ergonomic handle with clear imprint of the curve deflection, french size and sheath length
- Tip position indicator showing neutral or deflected curve



- Optimal SureSeal™ hemostatic valve technology to minimize air embolism and blood loss


- Dilator with depth markings for easy identification and placement
- Dilator cap with French size and guidewire indicator

TIP DEFLECTIONS



- 180° bi-directional deflection allows maximum REACH!
- Smooth tapered sheath tip to dilator transition for easy insertion during transseptal access
- Soft atraumatic tip designed to reduce potential for atrial trauma
- Radiopaque tip visibility enabling accurate tip positioning
- Reliable Flexsteer™ reinforced braid sheath technology provides superior control and kink resistance

ORDERING INFORMATION



Model # Set <i>(One per box)</i>	A1 Sheath Inner Diameter	A2 Usable Sheath Length	B1 Usable Dilator Length	B2 Overall Dilator Length	A3 Curve Bend	C Compatible Transseptal Needle
DCR0856717	● 8.5 F	67 cm	85 cm	88 cm	17 mm	89 cm
DCR0857317	● 8.5 F	73 cm	90 cm	93 cm	17 mm	98 cm
DCR0857717	● 8.5 F	77 cm	94 cm	97 cm	17 mm	98 cm
DCR0856722	● 8.5 F	67 cm	85 cm	88 cm	22 mm	89 cm
DCR0857322	● 8.5 F	73 cm	90 cm	93 cm	22 mm	98 cm
DCR0857722	● 8.5 F	77 cm	94 cm	97 cm	22 mm	98 cm
DCR0856750	● 8.5 F	67 cm	85 cm	88 cm	50 mm	89 cm
DCR0857350	● 8.5 F	73 cm	90 cm	93 cm	50 mm	98 cm
DCR0857750	● 8.5 F	77 cm	94 cm	97 cm	50 mm	98 cm
DCR1007317	● 10 F	73 cm	90 cm	93 cm	17 mm	98 cm
DCR1007717	● 10 F	77 cm	94 cm	97 cm	17 mm	98 cm
DCR1007322	● 10 F	73 cm	90 cm	93 cm	22 mm	98 cm
DCR1007722	● 10 F	77 cm	94 cm	97 cm	22 mm	98 cm
DCR1007350	● 10 F	73 cm	90 cm	93 cm	50 mm	98 cm
DCR1007750	● 10 F	77 cm	94 cm	97 cm	50 mm	98 cm
DCR1207122	● 12 F	71 cm	89 cm	93 cm	22 mm	98 cm
DCR1207522	● 12 F	75 cm	93 cm	97 cm	22 mm	98 cm
DCR1207150	● 12 F	71 cm	89 cm	93 cm	50 mm	98 cm
DCR1207550	● 12 F	75 cm	93 cm	97 cm	50 mm	98 cm

Guiding sheaths are supplied sterile. For single use only. Available bulk upon request.
Standard lengths are noted above. Other lengths available upon request.



P: 727.937.2511
F: 727.934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

