


Occlutech ASD Occluder

Compatibility chart & ordering information



Ø LA DISC: LEFT ATRIAL DISC DIAMETER
Ø WAIST: MAX. WAIST DIAMETER
Ø RA DISC: RIGHT ATRIAL DISC DIAMETER

Article No.	AS Defect Size (D) [mm]	Ø Waist [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	ODS* Size [F]	Pusher Art. No. (color)** Pistol Pusher / Flex Pusher
29ASD04	D≤4	4	11	9	7	55PP090 / 51FP060 ●
29ASD05	4<D≤5	5	14	11	7	55PP090 / 51FP060 ●
29ASD07	6<D≤7.5	7.5	18	14	7	55PP125 / 51FP100 ●
29ASD09	7.5<D≤9	9	20.5	16.5	7	55PP125 / 51FP100 ●
29ASD10	9<D≤10.5	10.5	22	18	7	55PP125 / 51FP100 ●
29ASD12	10.5<D≤12	12	27	23	9	55PP165 / 51FP120 ●
29ASD13	12<D≤13.5	13.5	28.5	24.5	9	55PP165 / 51FP120 ●
29ASD15	12<D≤15	15	30	26	9	55PP165 / 51FP120 ●
29ASD16	15<D≤16.5	16.5	31.5	27.5	9	55PP165 / 51FP120 ●
29ASD18	15<D≤18	18	33	29	9	55PP165 / 51FP120 ●
29ASD19	16.5<D≤19.5	19.5	34.5	30.5	10	55PP165 / 51FP120 ●
29ASD21	18<D≤21	21	36	32	11	55PP185 / 51FP150 ●
29ASD24	21<D≤24	24	39	35	11	55PP185 / 51FP150 ●
29ASD27	24<D≤27	27	42	38	12	55PP210 / 51FP160 ●
29ASD30	27<D≤30	30	45	41	12	55PP210 / 51FP160 ●
29ASD33	30<D≤33	33	48	43	12	55PP210 / 51FP160 ●
29ASD36	33<D≤36	36	52	46	12	55PP210 / 51FP160 ●
29ASD39	36<D≤39	39	54	49	12	55PP210 / 51FP160 ●
29ASD40	39<D≤40	40	55	50	12	55PP210 / 51FP160 ●

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

* The Occluders are compatible with Occlutech Delivery Set, e.g. ODS 45* (Art. nN. 51DSXXX). Recommended size. Sold separately.

**Color defines Pusher color



Occlutech Canada Inc.
2425 Matheson Blvd E, Suite 806
Mississauga ON L4W 5K4 Canada
order@occlutech.com, www.occlutech.com

