


COMPATIBILITY CHART

FIGULLA FLEX II PFO OCCLUDER



Ø LA: Left Atrial Disc Diameter,
Ø RA: Right Atrial Disc Diameter

Art. no.	PFO Size (D) [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	Flex-Pusher II (FP II) Art. no.	ODS* Size [F]
23PFO18DF	$D \leq 8$	16	18	51FP100 	7 F
23PFO25DF	$8 < D \leq 13$	23	25	51FP120 	9 F
23PFO30DF	$13 < D \leq 15$	27	30	51FP120 	9 F
23PFO35DF	$15 < D \leq 18$	31	35	51FP160 	11 F

* Figulla Flex II Occluder are compatible with Occlutech Delivery Set, e.g. ODS 45⁰ (Article no.: 51DSXXX)